

# **MYXOCOCCUS TOTAL RNA PREPARATION PROTOCOL**

*Adapted from existing protocols by Vinh Pham.*

*Last modified: February 12, 2002*

## **MATERIALS:**

**10% SDS**

**1M NaOAc pH 5.2**

**Phenol, Unbuffered/Water-saturated**

**24:1 Chloroform:Isoamyl Alcohol**

**3M NaOAc/1mM EDTA pH 8.0** (seal top with Parafilm to minimize evaporation)

**100% Ethanol** (store at -20°C)

**70% EtOH/10mM Tris pH 7.4/10mM NaCl** (seal top with Parafilm)

**DEPC-treated water**

**RQ1 DNase (1U/μl)** Promega cat# M6101

## **PROCEDURE:**

1. Grow cells to Klett 80-100. Pellet the desired amount of cells and resuspend in 500μl TPM. For developmental samples, spot 50-60 spots of 20μl each onto 150mm TPM agar plates and transfer fruiting bodies/spores/cells at the desired time point to a 1.5ml μfuge tube containing 500μl TPM.

*Note: The more cells, the greater the RNA yield. For volumes >2.0ml, pellet the cells in multiple 1.5ml μfuge tubes (1' at max speed) and resuspend in a single volume of 500μl TPM.*

2. Freeze at -80°C.

*Note: Ideally, cells should be quick frozen by immersing the μfuge tube in liquid nitrogen for ~1'.*

3. Pre-heat 10% SDS, 1M NaOAc pH 5.2, and phenol at 64°C.

4. Add 100μl SDS, mix, and shake every 10" for 1-2min. at 64°C.

*Note: Adding pre-heated SDS and pipetting to resuspend will thaw the sample.*

6. Add 110μl NaOAc.

7. Add 1ml phenol. Shake vigorously every 30"-1' for 4-6' at 64°C.

8. Spin on countertop microfuge for 10' at 14Krpm. Remove aqueous layer to new tube.

9. Repeat 1ml phenol extraction.

*Note: Perform additional extractions with phenol, as needed. Use phenol with a pH below 7.8 to minimize chromosomal DNA contamination. DNA enters the organic phase below this pH.*

10. Add 1ml RT chloroform/IAA. Spin for 10' at 14Krpm. Remove aqueous layer to new tube.

11. Add 130µl NaOAc-EDTA.
12. Add 1ml ice-cold 100% ethanol.

*Note: The sample can be stored at -20°C if desired.*

13. Spin for 25' at 10Krpm.
14. Wash pellet 3X with 500µl EtOH/Tris/NaCl.

*Note: Do not disrupt pellet between washes. Simply add the wash solution to the pellet, let sit for 1-5', then spin for 1' at max speed on a countertop µfuge.*

*Note: To prepare, autoclave*

<i>1M Tris, pH 7.6</i>	<i>1.76ml</i>
<i>NaCl</i>	<i>0.116g</i>
<i>dH<sub>2</sub>O</i>	<i>51.3ml</i>
<i>Adjust to pH 7.4</i>	

*Allow to cool.*

*Add 95% EtOH 147ml*

15. Resuspend in 100µl DEPC-treated water.

*Note: Add diethyl pyrocarbonate (DEPC) to a final concentration of 0.1% to distilled water, incubate at 37°C for 1hr, and inactivate by autoclaving.*

16. Add 5µl RQ1 DNase (RNase-free) & incubate at 37°C for 1hr to remove DNA contamination.

*Note: This step is essential if the RNA sample is to be used in quantitative experiments, e.g., RNA "slot blots" or real-time RT-PCR experiments.*

17. Incubate at 64°C for 5' to inactivate.