

Bradford Protein Assay (Bio-Rad)

Adapted from existing protocols by Vinh Pham

Last modified: June 22, 2003.

MATERIALS

Spectrophotometer

Immunoglobulin G, 1mg/ml

Resuspend lyophilized gamma globulin in dH₂O, store at -20°C

Disposable 1cm lightpath cuvettes

Polystyrene (340nm-750nm)

Methacrylate (275nm-750nm)

13 x 100 mm disposable borosilicate test tubes

Bio-Rad Protein Assay Dye (Bradford) Reagent Concentrate (5X)

PROCEDURE

1. Add the following volumes of the IgG protein standard to separate test tubes:

[IgG] _{final} in μ g/ml	Dilution Series	
	IgG (1mg/ml), μ l	dH ₂ O, μ l
0	0	100
2.5	2.5	97.5
5	5	95
7.5	7.5	92.5
10	10	90
12.5	12.5	87.5
15	15	85
17.5	17.5	82.5
20	20	80

Note: The linear range of the assay for IgG is 1.2 to 25 μ g/ml.

2. Prepare working solution of Bradford reagent. For each sample (including the IgG standards), prepare 0.9ml of working solution. The working solution is a mixture of Bradford reagent concentrate and water in a 2:7 ratio (e.g., 200 μ l Bradford reagent concentrate to 700 μ l water).

Note: Always prepare about twice as much working solution as is needed for each experimental sample in order to dilute samples with OD₅₉₅ readings outside the range of the IgG standard curve.

3. Prepare the experimental samples. Samples can be used as is, or they can be diluted in dH₂O. The final volume should be 100 μ l.

Note: Duplicate or triplicate samples should be prepared for more accurate [protein] measurements.

4. Transfer samples to separate test tubes.
5. Add 900 μ l Bradford reagent (working concentration) to each standard and experimental sample.

Note: The Coomassie Brilliant Blue G-250 dye shifts from 465nm to 595nm when binding to basic and aromatic amino acid residues, especially arginine (source: Bio-Rad Protein Assay manual).

6. Incubate at RT for at least 5 minutes.
7. Measure OD₅₉₅ absorption. If necessary, add excess Bradford reagent (plus 100 μ l dH₂O to make it a 1X working solution) prepared in step 2 to dilute samples to standard range.
8. Plot standard curve in a graphing utility (e.g., Microsoft Excel).
9. Plot a trendline (in Excel, go to Chart>Add Trendline...; click Options tab, check box for Display equation on chart). If not using a graphing utility, use $y=m(x-a)+b$ to calculate the equation for the trendline.
10. For each sample, the [protein] equals $(OD_{595}-b)/(1/m)(\text{dilution factor})$. The dilution factor is the final volume divided by the volume of protein sample (100 μ l or less).

Note: Adjust accordingly if excess working solution was added to dilute the sample. For example, if 5 μ l of protein sample was diluted with working solution to a final volume of 2ml (=2000 μ l), then the dilution factor is 2000/5, or 400.